

# BioComp Gradient Equipment Instruction

## Rotor SW41, Centrifuge Tube: SECTON 7030

### Laying gradients

- 1) Mark the halfway point of the tube using the marker block with the upper edge (short caps)
- 2) Light solution
  - a. Suck the light solution into a syringe (> 10 mL) with the provided cannula
  - b. Lay the light solution in the tube slightly above the line
- 3) Heavy solution
  - a. Suck the heavy solution into a syringe
  - b. Clean the outside of the cannula with paper towel
  - c. Push a drop of liquid out of the syringe and touch the cannula to paper towel to remove the excess solution
  - d. Push the cannula down the bottom of the tube along the wall
  - e. Laying the heavy solution and move the end of cannula upward to keep it 1 cm below the interface
  - f. When the interface reaches the line, move the cannula to the wall and pull the cannula out.
- 4) Cap the tube
  - a. Pivoting the cap into the tube so the side with the hole is the last to seal
  - b. Liquid should be visible inside the cap
- 5) Form the gradient
  - a. Turn on the system from the back
  - b. Select GMST
  - c. Use the UP and Down buttons to level the plate
  - d. Put tubes in the tube holder on the Gradient Master
  - e. Choose your gradient by selecting SW41 rotor; finding the gradient in the list; and hitting Run. Making a 10-50% gradient takes around 3 min
  - f. **No bubbles in the gradient!!! If there are bubbles, remake the gradient.**

*Users need to clean cannulas right after making gradients.*

### Loading samples

- 1) The gradient should be kept at 4 °C for at least 45 min before laying samples
- 2) Take the gradient cap off the centrifuge tube from the hole side first
- 3) After taking the cap off, the meniscus is about 3 mm below the opening. You should leave the **2 mm empty space after loading**. Therefore, you only can add 1 mm liquid, around 150 µL. To add more sample, you need to pipette out liquid from the top layer. For example, you want to **load the max amount 500 µL, you need to firstly pipette out 500-150=350 µL**.
- 4) Balance two tubes after pipetting out liquid.
- 5) After loading samples, confirm the weight difference, which should be less than 0.1 g.
- 6) Clean tubes outside before loading to a bucket

*Users need to clean gradient caps right after starting centrifugation.*

## Centrifugation

Follow your procedure to centrifuge samples

- 1) The balanced two tubes should be loaded opposite to each other on the rotor
- 2) All six buckets should be loaded to the rotor

## Fractionation

- 1) Fill up the rinse reservoir with MillQ water
- 2) Prime the pump
  - a. Push the toggle down and wait for water coming out of the waste tubing
  - b. Flip the toggle up
- 3) Turn the machine off and turn it back on
  - a. Hit the scan mode
- 4) Launch the software
  - a. Check the Core box
- 5) Turn on the fractionator
- 6) Screw the rinse adaptor to the piston
  - a. Push 5 ml water into system
- 7) Input Gradient information
- 8) Adjust the absorbance signal
  - a. Add water to the flow cell (Already done)
  - b. Adjust the On time to make the sample value between 800 K and 900 K
  - c. The On time should be change between run to run. If values change, the cell might not be cleaned well.
- 9) Zero OD (Only need Zero once a day)
  - a. Click the Start Zero OD button
- 10) Define the number of fractionating tubes
  - a. Select Meniscus Sensing
  - b. Change the collecting volume of each tube
- 11) Pre-Scan
  - a. Remove the rinse adaptor
  - b. Install the piston tips (need to be cleaned and dried between runs)**
- 12) Mount the tube in the holder
  - a. Push the mount cap on the tube and give a twist. The tube should be securely attached to the cap
  - b. Load the cap and the tube in to the mount.
  - c. Align the line on the cap with the screws on the mount and twist the cap clockwise to lock it
  - d. Slide the mount on the top of the light and twist clockwise
- 13) Name the run
- 14) Set up the fractionator
  - a. Load empty tubes on the tube rack
  - b. <30 tubes use the code 19 and load empty tubes to the middle two lanes,

- c. >30 tubes use the code 29 and load empty tubes to four lanes
- 15) Start scan from the software
- 16) Piston goes into the tube and run starts
- 17) Save Run

## **Cleaning**

- 1) Click the sterilization button
- 2) Rinse the piston
  - a. Take out the piston tip and connect the rinse adaptor
  - b. Wash the piston with 5-10 mL water by pushing the water through the adaptor
- 3) Rinse the system
  - a. Fill the rinse reservoir with MillQ water
  - b. Clean the system by holding Rinse and Air each 5 sec for 3 times
    - i. (Rinse 5 sec + Air 5 sec ) X 3
- 4) Sterilize the system
  - a. Empty water inside the rinse reservoir and fill it up with 15 mL 70% Ethanol
  - b. Rinse the system with 70% Ethanol
  - c. Hold air button for 5 sec
- 5) Air the system
  - a. Disconnect the rinse reservoir and connect an empty 60 mL syringe with the plunger pulling to the end
  - b. Air the pump
    - i. Push the toggle switch down
    - ii. Push air into the pump
    - iii. Retract the plunger and push air several times
  - c. Air the Triax
    - i. Push the toggle switch down + Retract the plunger
    - ii. Push the toggle switch up+ Push air
    - iii. Repeat step i and step ii three times